

SALTERSTAT™ CAPNOGRAPHY

Ordering Guide



SALTERSTAT™ REFLECTIVE SERIES


FOR MONITORS REQUIRING REFLECTIVE CONNECTORS



Oral-Nasal Facepiece




Nasal



Reflective Connector for Short-Term Use




Reflective Connector for Extended Use




Thread-Grip Connector



Nafion Tubing Used in Extended Use Codes



Filter



Soft Tubing

CO₂ SAMPLING WITH O₂ DELIVERY

ITEM	DURATION	INTERFACE	SIZE	O ₂ LENGTH	CO ₂ LENGTH	PK
4MSF1-0-6-25	Filtered for Short-Term Use	Oral-Nasal	Adult	2"	6'	25
4MSF1-0-6-100	Filtered for Short-Term Use	Oral-Nasal	Adult	2"	6'	100
4MSF1-0-13-25	Filtered for Short-Term Use	Oral-Nasal	Adult	2"	13'	25
4MSF1-7-6-25	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	6'	25
4MSF1-7-13-25	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	13'	25
4MSF1-7-6-100	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	6'	100
4MSF1-7-13-100	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	13'	100
4MSF2-PED-7-6-25	Filtered for Short-Term Use	Oral-Nasal	Pediatric	7'	6'	25
4MSF1-L-7-6-25	Filtered with Nafion for Extended Use	Oral-Nasal	Adult	7'	6'	25
4MSF1-L-7-13-25	Filtered with Nafion for Extended Use	Oral-Nasal	Adult	7'	13'	25
4MSF1-L-7-6-100	Filtered with Nafion for Extended Use	Oral-Nasal	Adult	7'	6'	100
4MSF2-L-PED-7-6-25	Filtered with Nafion for Extended Use	Oral-Nasal	Pediatric	7'	6'	25
4MSF2-L-PED-7-13-25	Filtered with Nafion for Extended Use	Oral-Nasal	Pediatric	7'	13'	25
4MSF3-7-6-25	Filtered for Short-Term Use	Nasal	Adult	7'	6'	25
4MSF3-7-6-100	Filtered for Short-Term Use	Nasal	Adult	7'	6'	100
4MSF3-7-13-25	Filtered for Short-Term Use	Nasal	Adult	7'	13'	25
4MSF3-7-13-100	Filtered for Short-Term Use	Nasal	Adult	7'	13'	100
4MSF3-PED-7-6-25	Filtered for Short-Term Use	Nasal	Pediatric	7'	6'	25
4MSF3-L-7-6-25	Filtered with Nafion for Extended Use	Nasal	Adult	7'	6'	25
4MSF3-L-PED-7-6-25	Filtered with Nafion for Extended Use	Nasal	Pediatric	7'	6'	25
4MSF3-L-INF-7-6-25	Filtered with Nafion for Extended Use	Nasal	Infant	7'	6'	25

CO₂ SAMPLING LINE

ITEM	DURATION	INTERFACE	SIZE	O ₂ LENGTH	CO ₂ LENGTH	PK
4MSF3-6-25	Filtered for Short-Term Use	Nasal	Adult	N/A	6'	25
4MSF3-INF-6-25	Filtered for Short-Term Use	Nasal	Infant	N/A	6'	25
4MSF3-L-6-25	Filtered with Nafion for Extended Use	Nasal	Adult	N/A	6'	25

SALTERSTAT™ REFLECTIVE SERIES

FOR MONITORS REQUIRING REFLECTIVE CONNECTORS



Reflective Series with Soft Tubing Headset



Sampling line with ETT Adapter

CO₂ SAMPLING WITH SOFT HEADSET TUBING

ITEM	DURATION	INTERFACE	SIZE	O ₂ LENGTH	CO ₂ LENGTH	PK
4MSF1SOFT-7-6-25	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	6	25
4MSF1SOFT-7-13-25	Filtered for Short-Term Use	Oral-Nasal	Adult	7'	13'	25
4MSF2SOFT-PED-7-6-25	Filtered for Short-Term Use	Oral-Nasal	Pediatric	7'	6	25
4MSF1SOFT-L-7-6-25	Filtered with Nafion for Extended Use	Oral-Nasal	Adult	7'	6	25
4MSF2SOFT-L-PED-7-6-25	Filtered with Nafion for Extended Use	Oral-Nasal	Pediatric	7'	6	25
4MSF3SOFT-7-6-25	Filtered with Nafion for Extended Use	Nasal	Adult	7'	6	25
4MSF3SOFT-7-13-25	Filtered with Nafion for Extended Use	Nasal	Adult	7'	13'	25
4MSF3SOFT-L-7-6-25	Filtered with Nafion for Extended Use	Nasal	Adult	7'	6	25

SAMPLING LINE WITH STANDARD ETT ADAPTER

ITEM	DURATION	SIZE	CO ₂ LENGTH	PK
4MSF5-6-25	Filtered for Short-Term Use	Adult/Pediatric	6'	25
4MSF5-6-100	Filtered for Short-Term Use	Adult/Pediatric	6	100
4MSF5-13-25	Filtered for Short-Term Use	Adult/Pediatric	13'	25
4MSF5-L-6-25	Filtered with Nafion for Extended Use	Adult/Pediatric	6'	25
4MSF5-L-6-100	Filtered with Nafion for Extended Use	Adult/Pediatric	6'	100
4MSF5-L-13-25	Filtered with Nafion for Extended Use	Adult/Pediatric	13'	25


CO₂ SAMPLING LINE WITH ONE REFLECTIVE, ONE LUER CONNECTOR

ITEM	DURATION	SIZE	CO ₂ LENGTH	PK
4MSF8-6-25	Filtered for Short-Term Use	Any	6'	25
4MSF8-13-25	Filtered for Short-Term Use	Any	13'	25
4MSF8-L-6-25	Filtered with Nafion for Extended Use	Any	6'	25

SALTERSTAT™ STANDARD SERIES

FOR MONITORS REQUIRING LUER CONNECTORS

CONNECTORS



Divided Cannula



Salter Eyes
*47 series, adult size only




CO₂/O₂ DIVIDED NASAL CANNULA WITH SOFT HEADSET TUBING

ITEM	SIZE	SALTER EYE	CO ₂ CONNECTOR	O ₂ CONNECTOR	CO ₂ LENGTH	O ₂ LENGTH	PK
47SOFT-7-0-25	Adult	Yes	Male Luer	Trumpet	2"	7'	25
47SOFT-7-7-25	Adult	Yes	Male Luer	Trumpet	7'	7'	25
47SOFT-F-7-0-25	Adult	Yes	Female Luer	Trumpet	2"	7'	25
47SOFT-F-7-7-25	Adult	Yes	Female Luer	Trumpet	7'	7'	25
47SOFT-FTG-10-10-25	Adult	Yes	Female Luer	Female Thread Grip	10'	10'	25
47SOFT-FTG-7-0-25	Adult	Yes	Female Luer	Female Thread Grip	2"	7'	25
47SOFT-FTG-7-7-25	Adult	Yes	Female Luer	Female Thread Grip	7'	7'	25
47SOFT-TG-14-14-25	Adult	Yes	Male Luer	Female Thread Grip	14'	14'	25
47SOFT-TG-7-0-25	Adult	Yes	Male Luer	Female Thread Grip	2"	7'	25
47SOFT-TG-7-7-25	Adult	Yes	Male Luer	Female Thread Grip	7'	7'	25
49SOFT-14-14-25	Adult	No	Male Luer	Trumpet	14'	14'	25
49SOFT-7-7-25	Adult	No	Male Luer	Trumpet	7'	7'	25
49SOFT-F-14-14-25	Adult	No	Female Luer	Trumpet	14'	14'	25
49SOFT-F-7-7-25	Adult	No	Female Luer	Trumpet	7'	7'	25
49SOFT-FTG-7-7-25	Adult	No	Female Luer	Female Thread Grip	7'	7'	25
49SOFT-TG-14-14-25	Adult	No	Male Luer	Female Thread Grip	14'	14'	25
49SOFT-TG-7-7-25	Adult	No	Male Luer	Female Thread Grip	7'	7'	25
49SOFT-FTG-10-10-25	Adult	No	Female Luer	Female Thread Grip	10'	10'	25
49SOFT-TG-10-10-25	Adult	No	Male Luer	Female Thread Grip	10'	10'	25
47SOFT-F-PED-14-14-25	Pediatric	No	Female Luer	Trumpet	14'	14'	25
47SOFT-FTG-PED-7-7-25	Pediatric	No	Female Luer	Female Thread Grip	7'	7'	25
47SOFT-PED-7-0-25	Pediatric	No	Male Luer	Trumpet	2"	7'	25
47SOFT-PED-7-7-25	Pediatric	No	Male Luer	Trumpet	7'	7'	25
47SOFT-TG-PED-7-7-25	Pediatric	No	Male Luer	Female Thread Grip	7'	7'	25
47SOFT-F-14-14-25	Adult	No	Female Luer	Trumpet	14'	14'	25

SALTERSTAT™ STANDARD SERIES

FOR MONITORS REQUIRING LUER CONNECTORS



CO₂/O₂ DIVIDED CANNULA WITH STANDARD HEADSET TUBING

ITEM	INTERFACE	SIZE	SALTER EYE	FILTER	CO ₂ CONNECTOR	O ₂ CONNECTOR	CO ₂ LENGTH	O ₂ LENGTH	PK
4700F-7-7-25	Nasal	Infant	Yes	No	Female Luer	Trumpet	7'	7'	25
4700FTG-7-7-25	Nasal	Infant	Yes	No	Female Luer	Female Thread Grip	7'	7'	25
4701F-7-0-25	Nasal	Infant	Yes	No	Female Luer	Trumpet	2"	7'	25
4706-7-0-25	Nasal	Adult	Yes	No	Male Luer	Trumpet	2"	7'	25
4706F-7-0-25	Nasal	Adult	Yes	No	Female Luer	Trumpet	2"	7'	25
4706FTG-7-0-25	Nasal	Adult	Yes	No	Female Luer	Female Thread Grip	2"	7'	25
4706TG-7-0-25	Nasal	Adult	Yes	No	Male Luer	Female Thread Grip	2"	7'	25
4707-10-10-25	Nasal	Adult	Yes	No	Male Luer	Trumpet	10'	10'	25
4707-13-13-25	Nasal	Adult	Yes	No	Male Luer	Trumpet	13'	13'	25
4707-7-7-25	Nasal	Adult	Yes	No	Male Luer	Trumpet	7'	7'	25
4707F-14-14-25	Nasal	Adult	Yes	No	Female Luer	Trumpet	14'	14'	25
4707F-7-7-25	Nasal	Adult	Yes	No	Female Luer	Trumpet	7'	7'	25
4707FTG-7-7-25	Nasal	Adult	Yes	No	Female Luer	Female Thread Grip	7'	7'	25
4707TG-14-14-25	Nasal	Adult	Yes	No	Male Luer	Female Thread Grip	14'	14'	25
4707TG-7-7-25	Nasal	Adult	Yes	No	Male Luer	Female Thread Grip	7'	7'	25
4950-10-10-25	Nasal	Adult	No	No	Male Luer	Trumpet	10'	10'	25
4950-14-14-25	Nasal	Adult	No	No	Male Luer	Trumpet	14'	14'	25
4950-7-7-25	Nasal	Adult	No	No	Male Luer	Trumpet	7'	7'	25
4950F-10-10-25	Nasal	Adult	No	No	Female Luer	Trumpet	10'	10'	25
4950F-14-14-25	Nasal	Adult	No	No	Female Luer	Trumpet	14'	14'	25
4950F-7-0-25	Nasal	Adult	No	No	Female Luer	Trumpet	2"	7'	25
4950F-7-7-25	Nasal	Adult	No	No	Female Luer	Trumpet	7'	7'	25
4950FTG-7-7-25	Nasal	Adult	No	No	Female Luer	Female Thread Grip	7'	7'	25
4950TG-14-14-25	Nasal	Adult	No	No	Male Luer	Female Thread Grip	14'	14'	25
4950TG-7-7-25	Nasal	Adult	No	No	Male Luer	Female Thread Grip	7'	7'	25
4700-7-7-25	Nasal	Infant	No	No	Male Luer	Trumpet	7'	7'	25
4701-7-0-25	Nasal	Infant	No	No	Male Luer	Trumpet	2"	7'	25
4702-7-0-25	Nasal	Pediatric	No	No	Male Luer	Trumpet	2"	7'	25
4702F-7-0-25	Nasal	Pediatric	No	No	Female Luer	Trumpet	2"	7'	25
4703-7-7-25	Nasal	Pediatric	No	No	Male Luer	Trumpet	7'	7'	25
4703F-7-7-25	Nasal	Pediatric	No	No	Female Luer	Trumpet	7'	7'	25
4703FTG-7-7-25	Nasal	Pediatric	No	No	Female Luer	Female Thread Grip	7'	7'	25
4703TG-7-7-25	Nasal	Pediatric	No	No	Male Luer	Female Thread Grip	7'	7'	25
4704-7-7-25	Nasal	Infant	No	No	Male Luer	Trumpet	7'	7'	25
4704F-7-7-25	Nasal	Infant	No	No	Female Luer	Trumpet	7'	7'	25
4705-7-0-25	Nasal	Infant	No	No	Male Luer	Trumpet	2"	7'	25
4705F-7-0-25	Nasal	Infant	No	No	Female Luer	Trumpet	2"	7'	25
4793-7-7-25	Oral-Trac Trunk	Pediatric	No	Yes	Male Luer	Trumpet	7'	7'	25
4793F-7-7-25	Oral-Trac Trunk	Pediatric	No	Yes	Female Luer	Trumpet	7'	7'	25

SALTERSTAT™ STANDARD SERIES

FOR MONITORS REQUIRING LUER CONNECTORS



Piggyback



CO₂ Nasal Sampling



Oral-Trach®

PIGGYBACK CO₂/O₂ NASAL SAMPLING CANNULAS

ITEM	INTERFACE	SIZE	SALTER EYE	O ₂ CONNECTOR	CO ₂ CONNECTOR	CO ₂ LENGTH	O ₂ LENGTH	PK
4002-7-7-25	Nasal	Adult	No	Trumpet	Male Luer	7'	7'	25
4002-14-14-25	Nasal	Adult	No	Trumpet	Male Luer	14'	14'	25
4002F-7-7-25	Nasal	Adult	No	Trumpet	Female Luer	7'	7'	25
4102-7-7-25	Nasal	Pediatric	No	Trumpet	Male Luer	7'	7'	25
4102F-7-7-25	Nasal	Pediatric	No	Trumpet	Female Luer	7'	7'	25

CO₂ NASAL SAMPLING CANNULAS (NO O₂ DELIVERY)

ITEM	INTERFACE	SIZE	CO ₂ CONNECTOR	O ₂ LENGTH	PK
4000-7-25	Nasal	Adult	Male Luer	7'	25
4000F-7-25	Nasal	Adult	Female Luer	7'	25
4100-7-25	Nasal	Pediatric	Male Luer	7'	25
4007-7-25	Nasal	Adult	Male Luer	7'	25
4100F-7-25	Nasal	Pediatric	Female Luer	7'	25
4001-7-25	Oral-Nasal	Adult	Male Luer	7'	25
4101-7-25	Oral-Nasal	Pediatric	Male Luer	7'	25
4200-7-27	Nasal	Infant	Male-Luer	7'	25

ORAL-TRAC® ORAL-NASAL CANNULA

ITEM	INTERFACE	SIZE	O ₂ CONNECTOR	CO ₂ CONNECTOR	CO ₂ LENGTH	O ₂ LENGTH	PK
4797TG-7-7-25	Oral-Trac Trunk	Adult	Female Thread Grip	Male Luer	7'	7'	25
4797-7-7-25	Oral-Trac Trunk	Adult	Trumpet	Male Luer	7'	7'	25
4797F-7-7-25	Oral-Trac Trunk	Adult	Trumpet	Female Luer	7'	7'	25
4796F-7-0-25	Oral-Trac Trunk	Adult	Trumpet	Female Luer	2'	7'	25

GAS SAMPLING LINES

ITEM	DESCRIPTION	MOISTURE MANAGEMENT	CONNECTOR	TUBING I.D.	LENGTH	PK
0512	Standard	No	2 Male Luer	0.05"	10'	50
0513	Standard	No	1 Male, 1 Female Luer	0.05"	10'	50
0516	Co-extruded PVC/Polyethylene	No	2 Male Luer	0.047"	10'	50
0533	Standard	No	2 Male Luer	0.05"	15'	50
0535	Standard with Swivel Elbow	No	2 Male Luer	0.05"	10'	50
ZML-9333-8-50	Standard	12" Nafion	2 Male Luer	0.05"	8'	50
ZML-9336-8-50	Standard	12" Nafion	2 Male Luer	0.05"	8'	50
ZML-9530B-10-500	Co-Extruded PVC/Polyethylene	No	2 Male Luer	0.05"	10'	500
ZML-9540B-10-500	Co-Extruded PVC/Polyethylene	No	2 Male Luer	0.047"	10'	500
ZML-9800-1-50	Standard	12" Nafion	1 Male, 1 Female Luer	0.05"	8'	50
ZML-9801-1-25	Standard	6" Nafion	1 Male, 1 Female Luer	0.05"	8'	25
ZML-9803-4-25	Standard	48" Nafion	1 Male, 1 Female Luer	0.05"	8'	25

CAPNOVUE® MASKS



CapnoVue M1 Mask



CapnoVue Scope Mask

CAPNOVUE® CAPNOGRAPHY MASKS

ITEM	MASK STYLE	CONNECTOR	SIZE	LENGTH	PK
4727-7-0-50	M1 Mask without Oral Access	Trumpet	Adult	7'	50
4727-14-0-50	M1 Mask without Oral Access	Trumpet	Adult	14'	50
4723-7-0-50	M1 Mask without Oral Access	Trumpet	Pediatric	7'	50
4727TG-7-0-50	M1 Mask without Oral Access	Female Thread Grip	Adult	7'	50
4727TG-14-0-50	M1 Mask without Oral Access	Female Thread Grip	Adult	14'	50
4723TG-7-0-50	M1 Mask without Oral Access	Female Thread Grip	Pediatric	7'	50
4721TG-7-0-50	M1 Mask without Oral Access	Female Thread Grip	Infant	7'	50
4737-7-0-50	Scope Mask with Oral Access	Trumpet	Adult	7'	50
4737-14-0-50	Scope Mask with Oral Access	Trumpet	Adult	14'	50
4733-7-0-50	Scope Mask with Oral Access	Trumpet	Pediatric	7'	50
4737TG-7-0-50	Scope Mask with Oral Access	Female Thread Grip	Adult	7'	50
4733TG-7-0-50	Scope Mask with Oral Access	Female Thread Grip	Pediatric	7'	50
4745-14-15-50	Scope Mask with Oral Access	Luer	Adult	14'	50
4745-7-7-50	Scope Mask with Oral Access	Luer	Adult	7'	50
47MSF8-14-13-50	Scope Mask with Oral Access	Reflective	Adult	14'	50

All CapnoVue® Masks packaged 50 per case. CO₂ Sampling line is not included. Nafion® is a registered trademark of Dupont, CapnoVue® is a registered trademark of Salter Labs®, Luer-Lok® and Luer-Slip® are registered trademarks of Becton Dickinson



GDE-20002-C1 Rev1

